INOVO
SENSE

SHAPE Plus

INSTRUCTIONS FOR USE

1. Description

• **INOVOSENSE SHAPE Plus** is a sterile, non-pyrogenic, viscoelastic, colourless, transparent cross-linked sodium hyaluronate gel of non-animal origin in a physiological phosphate buffer. **INOVOSENSE SHAPE Plus** contains 0.3% lidocaine hydrochloride.

2. Presentation

• **INOVOSENSE SHAPE Plus** is presented in a single use pre-filled glass syringe sterilized by moist heat. Each box contains one instruction leaflet, two syringes, two traceability labels and six sterile needles, two micro cannulas for single use only. The dimensions of needles and micro cannulas are stated on the external box.

3. Composition

Cross-linked sodium hyaluronate: 25 mg

Lidocaine hydrochloride: 3.0 mg

Phosphate buffer pH 7 q.s.: 1ml

4. Intended use/Indications

4.1 Intended use

• **INOVOSENSE SHAPE Plus** is an injectable biodegradable implant intended for restoration of facial volume and also correcting signs of facial fat loss (lipoatrophy) in patients with HIV.

• **INOVOSENSE SHAPE Plus** is also intended for restoration of facial volume.

4.2 Indications

• **INOVOSENSE SHAPE Plus** is indicated for the restoration or augmentation of facial volume by injection into the deep dermis or

subcutaneously.

- The presence of lidocaine aims to reduce local pain associated with the injection of the gel and to improve patient comfort.

5. Posology and administration method

- **INOVOSENSE SHAPE Plus** is designed to be injected into deep subcutaneous layers or above the periosteum by a legally authorized practitioner who have the appropriate training, experience and who are knowledgeable about the anatomy at and around the site of injection in order to minimize the risk of potential complications.

- Inject **INOVOSENSE SHAPE Plus** slowly and not too fast to apply the least amount of pressure necessary, according to the appropriate injection technique using the provided needles.

- General recommended injection techniques are for example: linear or serial threading, fanning, cross-hatching or serial puncture. The quantity of product to be injected depends on the area to be corrected.

- The risk of an intravascular injection can be reduced by different strategies, including aspiration prior to injection, utilizing lower volumes and serial injections in high-risk areas, treating one side at a time, pinching/tenting the skin to provide more space superficial to the branches of the main arteries, and manual occlusion of the origin of the supratrochlear vessels with the

non-dominant finger. Blunt cannulas may reduce, but not eliminate the risk.

- **INOVOSENSE SHAPE Plus** can be used for all Fitzpatrick skin types.
- **INOVOSENSE SHAPE Plus** must be injected under appropriate aseptic conditions into healthy, non-inflamed skin. Before injection, thoroughly disinfect the area to be treated.
- To ensure optimal use of **INOVOSENSE SHAPE Plus**, it is recommended to assemble the needle according to the diagrams below. Improper assembly may lead to a separation of the needle and syringe and/or leakage of the gel at the Luer-lock connection during injection.
- If the needle becomes obstructed and the injection pressure is too high, stop the injection and change the needle.
- The use of the enclosed needle is recommended, as a smaller needle diameter would require a greater force to inject the implant.
- The quantity of the gel to be injected depends on the area to be treated and the correction to be achieved. Do not over-correct.
- Gently massage the treated area after the injection to distribute the product uniformly.

6. Contra-indications

- **INOVOSENSE SHAPE Plus** is contra-indicated:
 - In case of known hypersensitivity to one of the product's components, especially to sodium hyaluronate, BDDE, lidocaine hydrochloride and amide-type local anesthetics,

- In pregnant and breast-feeding women,
 - In patients under 18 years old,
 - In patients presenting a general infection,
 - In patients presenting an active auto-immune disease.
- Do not inject **INOVOSENSE SHAPE Plus** into blood vessels.
 - Do not inject **INOVOSENSE SHAPE Plus** into skin areas presenting active cutaneous inflammation or infection due to e.g. immunological, allergic, bacterial, fungal or viral causes.
 - Do not inject **INOVOSENSE SHAPE Plus** into an area previously treated with a permanent dermal filler.
 - Do not inject **INOVOSENSE SHAPE Plus** into correction of superficial wrinkles and fine lines (injection of the product into the superficial dermis).

7. Precautions for use

- Healthcare practitioners are encouraged to discuss all potential risks of soft tissue injection with their patients prior to treatment and ensure that patients are aware of signs and symptoms of potential complications.
- In the absence of available clinical data on tolerance of the injection of **INOVOSENSE SHAPE Plus** in patients presenting a history of severe multiple allergies or anaphylactic shock, the practitioner must decide whether to inject **INOVOSENSE SHAPE Plus** on a case-by-case basis depending on the nature of the disease as well as the

associated treatment as it may worsen the existing patient health condition. It is recommended to propose a prior double test to these patients and not to inject if the disease is evolving. It is also recommended to carefully monitor these patients after injection.

- It is recommended not to inject **INOVOSENSE SHAPE Plus** in patients with a history of streptococcal diseases and in patients pre-disposed to hypertrophic scars or keloids.

- **INOVOSENSE SHAPE Plus** injected in the temples area may be associated with an increased risk for intravascular complications and the consequences of local vascular occlusion, embolization, vision impairment, blindness, ischemia, necrosis or infarction.

- **INOVOSENSE SHAPE Plus** can be used in combination with other **INOVOSENSE** products during the same session but in different facial areas. Instructions for use of each product should be followed.

- **INOVOSENSE SHAPE Plus** can be used in combination treatments such as with botulinum toxin and/or calcium hydroxylapatite (Radiesse®) only if injected in different facial areas. Practitioners should be experienced and patients appropriately selected as benefits but also adverse events can be cumulative and causality of adverse events could become difficult to determine. Instructions for use, depth of injection and appropriate recommendation of each product should be followed.

- No clinical data are available on the injection of **INOVOSENSE SHAPE Plus** into an area already treated with other aesthetic products or procedures.
- **INOVOSENSE SHAPE Plus** must not be used in association with other aesthetic techniques such as peeling, dermabrasion or any type of laser treatment before complete healing of the last treatment. In any case, even if the healing occurs earlier, **INOVOSENSE SHAPE Plus** must not be used earlier than 2 weeks after the last treatment.
- Patients using anti-coagulation, anti-platelet, or thrombolytic medications (e.g. warfarin), anti-inflammatory drugs (oral/injectable corticosteroids or non-steroidal anti-inflammatory drugs (NSAIDs; e.g. aspirin, ibuprofen) or other substances known to increase coagulation time (vitamins or herbal supplements, e.g. Vitamin E, garlic, Ginkgo biloba and St. John's Wort), from 10 days pre- to 3 days post-injection may have increased reactions of hematomas, nodules or bleeding at the injection site.
- Injection of **INOVOSENSE SHAPE Plus** into patients with a history of previous herpetic eruption may be associated with reactivation of the herpes (HHV related diseases, e.g. pityriasis rosea).
- In cases of patients suffering from epilepsy, impaired cardiac conditions, severely impaired hepatic function or severe renal dysfunction or porphyria, the practitioner must de-

cide whether to inject **INOVOSENSE SHAPE Plus** on a case-by-case basis depending on the nature of the disease as well as the associated treatment.

- Practitioners and athletes should consider that lidocaine may produce positive results to anti-doping tests.
- It should be noted that the presence of lidocaine may cause local redness or hypersensitivity.
- For normal healthy adults it is recommended that the maximum total dose of lidocaine HCl (without epinephrine) does not exceed 300 mg (or 4.5 mg/kg) per session. Overdosage of lidocaine HCl usually results in sign of the central nervous system or cardiovascular toxicity.
- When using concurrently (topical administration...), the total administered dose of lidocaine should be considered.
- Care should be taken for patients with congenital methemoglobinemia, with glucose-6-phosphate dehydrogenase deficiencies and patients who are receiving concomitant treatment with methemoglobin- inducing agents.
- Check the integrity of the inner packaging and the expiry date for both the syringe and the needle prior to use. Do not use these products if the expiry date has lapsed or if the inner packaging has been opened or damaged.
- Do not transfer **INOVOSENSE SHAPE Plus** into another container and do not add other substances to the product.

- Only the gel is sterile, but not the outside of the syringe.
- Discard the syringe, the remaining product and the needles in the appropriate container after use.
- Do not re-sterilize and do not reuse due to the associated risks including infection.
- The patient must avoid applying makeup (including skin care products) for at least 12 hours after treatment as well as avoid saunas, peeling, Turkish baths and prolonged exposure to the sun, UV rays, extreme heat and cold for 2 weeks after the treatment. Patients should also avoid putting pressure on and/or handling the treated area and should avoid strenuous physical activity following treatment.
- The patient must avoid drinking alcohol for 24 hours before and after treatment. Alcohol may cause the blood vessels to dilate and cause more bruising.

8. Warnings

- Sodium hyaluronate precipitates in the presence of quaternary ammonium salts (such as benzalkonium chloride). It is therefore recommended that **INOVOSENSE SHAPE Plus** does not come into contact with such substances.
- Rare but serious adverse events associated with the intravascular injection of soft tissue fillers in the face have been reported and include temporary or permanent vascular complication,

vision impairment, blindness, cerebral ischemia or cerebral hemorrhage, leading to stroke, skin necrosis, and damage to underlying facial structures. Practitioners should immediately stop the injection if a patient exhibits any of the following symptoms, including changes in vision, signs of a stroke, blanching of the skin or unusual pain during or shortly after the procedure.

- Patients should receive prompt medical attention and possibly evaluation by an appropriate healthcare practitioner specialist should an intravascular injection occur.

9. Side effects and adverse events

- Patients must be informed by the practitioner about possible side effects and adverse events before treatment.

• Side effects

- Injection site reactions may occur following injection into the skin but disappear spontaneously within a few days. This includes swelling, nodule or lump/bump, bruising, hematoma, ecchymosis, induration, erythema/redness, tenderness, pain, discoloration, pruritus/itching, tingling, paraesthesia, numbness, hypoesthesia, scabbing, needle mark and discomfort or irritation.
- These injection site reactions are generally of mild or moderate intensity. A transient bleeding may also occur at the injection site and usually stops spontaneously as soon as the injection is finished.

• Adverse events

• In occasional cases, one or more of the following may occur either immediately or as a delayed reaction: acne cystic, milia, skin dryness (rough facial skin, exfoliation), injection site erosion, inflammation, shivering, fatigue, lymphatic system disorder, rash, burning sensation, injection site sore/warmth/fever, pruritus/itching, urticaria, hematoma, telangiectasia, ecchymosis, edema (including lymph edema), headache/cephalgia, tumefaction, tension, swelling (including persistent swelling), hyper- or hypo-pigmentation, angioedema, induration, blister vesicle, papule, lump/bump (visible and/or palpable material) or nodule (including inflammatory nodules), mass, granuloma (including inflammatory signs and foreign body reactions), necrosis, ischemia, vascular occlusion, embolization, infarction, Tyndall effect (including translucent chords), hypersensitivity, allergic reactions (including asthma attack, Quincke's edema, anaphylactic shock or throat tightening) to one of the product's components (e.g. hyaluronic acid, BDDE), oral and dental disorders, nervous system impairment, impairment of the otorhinolaryngological system (e.g. nasal congestion, oropharyngeal pain, dysgeusia, rhinorrhea, epistaxis, sinusitis, transient hearing loss), mastication pain, parotid gland enlargement, muscle twitching, muscle injury, nausea, vomiting, circulatory collapse, presyncope, peripheral venous disease, hot flush, anxiety caused by trypanophobia, patient dissatisfaction and disappointment (due to lack of or reduced performance, decreased

firmness/response, undesirable aesthetic effect), injection site discharge, device migration, product distribution issue (e.g. product accumulation, injection site indentation, superficial vein prominence, overcorrection or cranial nerve disorder (e.g. cranial nerve paralysis, facial paralysis, trigeminal neuralgia).

- Literature reported rare cases of the following adverse events with hyaluronic acid products such as infection (e.g. erysipelas, phlegmon, cellulitis, including open or draining wounds and (dental) abscess, impetigo, pustules), chronic infection (including biofilm formation), scarring, persistent skin discoloration, sensory dysfunction, non-thrombotic lung embolism as well as sarcoid granuloma formation in subjects with hepatitis C and interferon treatment, cerebral injuries (e.g. intracranial penetration, subarachnoid hemorrhage), strabismus, ophthalmoplegia, iris adhesions, cataract, conjunctival hemorrhage, eyelid ptosis and lacrimation.

- The risk of granuloma, ischemia, necrosis and vascular occlusion is higher with deep injections and high volumes.

- Patients should be instructed to report any side effects which last for more than one week and any adverse event as soon as it occurs to his/her practitioner, especially if patient has changes in his/her vision, signs of a stroke (including sudden difficulty speaking, numbness or weakness in his/her face, arms, or legs, difficulty

walking, face drooping, severe headache, dizziness, or confusion), white appearance of the skin, or unusual pain during or shortly after treatment.

- The practitioner may then refer the patient to the appropriate treatment.

- Patients with lighter skin types are more likely to develop injection-related adverse events. However, patients with skin of color are more likely to develop post-inflammatory hyperpigmentation and/or hypertrophic scar/keloid formation following injection procedures. Patients with specific ethnic characteristics, e.g. Asian population, should be informed of a higher risk of tissue reactions, e.g. itching, swelling, erythema, inflammation.

10. Assembly of needle to syringe

- For optimal use of **INOVOSENSE SHAPE Plus**, it is important that the needle is properly connected to the syringe. See diagrams 1, 2, 3 and 4.

1. **Firmly hold** the glass cylinder of **the syringe and the Luer-lock** adaptor between the thumb and the index fingers.





2. Grasp the protective cap with the other hand and break the seal.

3. **Push & Twist** the needle on the syringe **until a resistance is felt**. Do not over-tight. Overtightening of the needle may lead to the Luer-lock moving and dislodging from the syringe.

4. Keep holding the Luer-lock and remove the sheath from the needle.

11. Storage

Store between 2°C and 25°C. Protect from light and freezing. Avoid mechanical shocks.

Marketing Authorization Holder:

MEDISA ARA GOSTAR IRANIAN

Number 21, Jooybar Alley Shad St,

Molla Sadra St, Tehran, Iran

Manufacturer:

ARYA TAJHIZ IBN SINA

Number 21, 20 St, Kaveh Industrial City, Saveh, Iran

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