INOVO
SENSE

SHINE Plus

INSTRUCTIONS FOR USE

1. Description

. **INOVOSENSE SHINE Plus** is a sterile, non-pyrogenic, viscoelastic, colourless, transparent cross-linked sodium hyaluronate gel of non-animal origin in a physiological phosphate buffer. **INOVOSENSE SHINE Plus** contains 0.3% lidocaine hydrochloride.

2. Presentation

. **INOVOSENSE SHINE Plus** is presented in a single use pre-filled glass syringe sterilized by moist heat. Each box contains one instruction leaflet, two syringes, two traceability labels and four sterile needles for single use only. The dimensions of needles are stated on the external box.

3. Composition

Cross-linked sodium hyaluronate: 10 mg

Lidocaine hydrochloride: 3.0 mg

Phosphate buffer pH 7 q.s.: 1ml

4. Intended use/Indications

4.1 Intended use

. **INOVOSENSE SHINE Plus** is an injectable biodegradable implant intended for filling of depressions due to facial atrophic scars and for skin re-hydration and rejuvenation.

4.2 Indications

. **INOVOSENSE SHINE Plus** is indicated for filling of depressions due to facial atrophic scars and for skin re-hydration and rejuvenation. **INOVOSENSE SHINE Plus** enables to increase skin firmness and elasticity and to improve

complexion and radiance of the face, neck, décolletage and hands.

. The presence of lidocaine aims to reduce local pain associated with the injection of the gel and to improve patient comfort.

5. Posology and administration method

. **INOVOSENSE SHINE Plus** is designed to be injected into the superficial dermis or dermo-epidermal junction by a legally authorized practitioner who have the appropriate training, experience, and who are knowledgeable about the anatomy at and around the site of injection in order to minimize the risk of potential complications.

. Facial atrophic scars are treated through one treatment session in the superficial part of the dermis. Skin rehydration and rejuvenation can be achieved with a treatment plan of 3 sessions of multipuncture or micro-papular injections 4 weeks apart into the superficial dermis or the dermo-epidermal junction.

. **INOVOSENSE SHINE Plus** may be injected using the multipuncture or the micro-papular techniques. It is advisable to use the multipuncture technique over the whole face performing a double series of parallel injections, the 2 series being perpendicular to each other. The micro-papules obtained by repeated injections of small volumes can be considered on the most dehydrated areas. If this latter technique is used, it is recommended that the patient be informed that the micro-pap-

ules may be visible for 2 to 3 days.

- **INOVOSENSE SHINE Plus** must be injected into healthy, non-inflamed and previously appropriately disinfected skin.

- It is recommended to use the supplied needles.

To ensure optimal use of **INOVOSENSE SHINE Plus** it is recommended to assemble the needle according to the diagrams below.

- If the needle becomes obstructed and the injection pressure is too high, stop the injection and change the needle.

- **INOVOSENSE SHINE Plus** should be injected slowly. The quantity of the gel to be injected depends on the area to be treated. Gently massage the treated area after the injection to distribute the product uniformly.

6. Contra-indications

- **INOVOSENSE SHINE Plus** is contra-indicated:

- In case of known hypersensitivity to one of the product's components, especially to sodium hyaluronate, BDDE, lidocaine hydrochloride and amide-type local anesthetics,

- In pregnant and breast-feeding women,

- In patients under 18 years old,

- In patients presenting a general infection

- In patients presenting an active auto-immune disease.

- Do not inject **INOVOSENSE SHINE Plus** into blood vessels.

- Do not inject **INOVOSENSE SHINE Plus** into skin areas presenting active cutaneous inflammation

or infection due to e.g. immunological, allergic, bacterial, fungal or viral causes.

- Do not inject **INOVOSENSE SHINE Plus** into an area previously treated with a permanent dermal filler.

7. Precautions for use

- Healthcare practitioners are encouraged to discuss all potential risks of soft tissue injection with their patients prior to treatment and ensure that patients are aware of signs and symptoms of potential complications.

- In the absence of available clinical data on tolerance of the injection of **INOVOSENSE SHINE Plus** in patients presenting a history of severe multiple allergies or anaphylactic shock, the practitioner must decide whether to inject **INOVOSENSE SHINE Plus** on a case-by-case basis depending on the nature of the disease as well as the associated treatment as it may worsen the existing patient health condition. It is recommended to propose a prior double test to these patients and not to inject if the disease is evolving. It is also recommended to carefully monitor these patients after injection.

- It is recommended not to inject **INOVOSENSE SHINE Plus** in patients with a history of streptococcal diseases and in patients pre-disposed to hypertrophic scars or keloids.

- **INOVOSENSE SHINE Plus** can be injected in the periorbital region including tear trough only by trained experienced practitioners with a deep knowledge of the anatomy. Injection of

soft tissue gels into this area may be associated with an increased frequency and severity of side effects.

- **INOVOSENSE SHINE Plus** can be used in combination with other **INOVOSENSE** products during the same session but in different facial areas. Instructions for use of each product should be followed.

- **INOVOSENSE SHINE Plus** can be used in combination treatments such as with botulinum toxin and/or calcium hydroxylapatite (Radiesse®) only if injected in different facial areas. Practitioners should be experienced and patients appropriately selected as benefits but also adverse events can be cumulative and causality of adverse events could become difficult to determine. Instructions for use, depth of injection and appropriate recommendation of each product should be followed. No clinical data are available on the injection of **INOVOSENSE SHINE Plus** into an area already treated with other aesthetic products or procedures.

- **INOVOSENSE SHINE Plus** must not be used in association with other aesthetic techniques such as peeling, dermabrasion or any type of laser treatment before complete healing of the last treatment. In any case, even if the healing occurs earlier, **INOVOSENSE SHINE Plus** must not be used earlier than 2 weeks after the last treatment.

- **INOVOSENSE SHINE Plus** can be used in addition with laser treatments that should be conducted 1 or 2 months prior to hyaluronic acid injection.

Check the integrity of the inner packaging and the expiry date for both the syringe and the needle prior to use. Do not use these products if the expiry date has lapsed or if the inner packaging has been opened or damaged.

- Patients using anti-coagulation, anti-platelet, or thrombolytic medications (e.g. warfarin), anti-inflammatory drugs (oral/injectable) corticosteroids or non-steroidal anti-inflammatory drugs (NSAIDs; e.g. aspirin, ibuprofen), or other substances known to increase coagulation time (vitamins or herbal supplements, e.g. Vitamin E, garlic, Ginkgo biloba and St. John's Wort), from 10 days pre- to 3 days post-injection may have increased reactions of hematomas, nodules or bleeding at the injection site.

- Injection of **INOVOSENSE SHINE Plus** into patients with a history of previous herpetic eruption may be associated with reactivation of the herpes (HHV related diseases, e.g. pityriasis rosea).

- In cases of patients suffering from epilepsy, impaired cardiac conditions, severely impaired hepatic function or severe renal dysfunction or porphyria, the practitioner must decide whether to inject **INOVOSENSE SHINE Plus** on a case-by-case basis depending on the nature of the disease as well as the associated treatment.

Practitioners and athletes should consider that lidocaine may produce positive results to anti-doping tests.

- It should be noted that the presence of lidocaine

may cause local redness or hypersensitivity.

For normal healthy adults it is recommended that the maximum total dose of lidocaine HCl (without epinephrine) does not exceed 300 mg (or 4.5 mg/kg) per session. Overdosage of lidocaine HCl usually results in sign of the central nervous system or cardiovascular toxicity.

- When using concurrently (topical administration...), the total administered dose of lidocaine should be considered.

- Care should be taken for patients with congenital methemoglobinemia, with glucose-6-phosphate dehydrogenase deficiencies and patients who are receiving concomitant treatment with methemoglobin-inducing agents.

- Check the integrity of the inner packaging and the expiry date for both the syringe and the needle prior to use. Do not use these products if the expiry date has lapsed or if the inner packaging has been opened or damaged.

- Do not transfer **INOVOSENSE SHINE Plus** into another container and do not add other substances to the product.

- Only the gel is sterile, but not the outside of the syringe.

- Discard the syringe, the remaining product and the needles in the appropriate container after use.

- Do not re-sterilize and do not reuse due to the associated risks including infection.

- The patient must avoid applying makeup (including skin care products) for at least 12 hours after

treatment as well as avoid saunas, peeling, Turkish baths and prolonged exposure to the sun, UV rays, extreme heat and cold for 2 weeks after the treatment. Patients should also avoid putting pressure on and/or handling the treated area and should avoid strenuous physical activity following treatment.

. Patients should also avoid putting pressure on and/or handling the treated area. When the micro-papular technique has been used, the patient can massage the papules at the condition that the hands and treated area are disinfected.

8. Warnings

· Sodium hyaluronate precipitates in the presence of quaternary ammonium salts (such as benzalkonium chloride). It is therefore recommended that **INOVOSENSE SHINE Plus** does not come into contact with such substances.

· Rare but serious adverse events associated with the intravascular injection of soft tissue fillers in the face have been reported and include temporary or permanent vascular complication, vision impairment, blindness, cerebral ischemia or cerebral hemorrhage, leading to stroke, skin necrosis, and damage to underlying facial structures. Practitioners should immediately stop the injection if a patient exhibits any of the following symptoms, including changes in vision, signs of a stroke, blanching of the skin or unusual pain during or shortly after the procedure. Patients should receive prompt medical attention and possibly evaluation

by an appropriate healthcare practitioner specialist should an intravascular injection occur.

9. Side effects and adverse events

. Patients must be informed by the practitioner about possible side effects and adverse events before treatment.

• Side effects

. Injection site reactions may occur following injection into the skin but disappear spontaneously within a few days. This includes swelling, nodule or lump/bump, bruising, hematoma, ecchymosis, induration, erythema/redness, tenderness, pain, discoloration, pruritus/itching, tingling, paraesthesia, numbness, hypoesthesia, scabbing, needle mark and discomfort or irritation. These injection site reactions are generally of mild or moderate intensity. A transient bleeding may also occur at the injection site and usually stops spontaneously as soon as the injection is finished.

• Adverse events

In occasional cases, one or more of the following may occur either immediately or as a delayed reaction: inflammation, rash, burning sensation, pruritus/itching, hematoma, edema, swelling, hyper- or hypo-pigmentation, Tyndall effect, blister, induration, lump/ bump or nodule.

. Literature reported rare cases of the following adverse events with hyaluronic acid products such as infection and abscess, scarring, hypersensitivity or allergy to one of the product's components, granuloma, ischemia or necrosis. The risk is higher

with deep injections and high volumes.

. Patients should be instructed to report any side effects which last for more than one week and any adverse event as soon as it occurs to his/her practitioner, especially if patient has changes in his/her vision, signs of a stroke (including sudden difficulty speaking, numbness or weakness in his/ her face, arms, or legs, difficulty walking, face drooping, severe headache, dizziness, or confusion), white appearance of the skin, or unusual pain during or shortly after treatment. The practitioner may then refer the patient to the appropriate treatment.

. Patients with lighter skin types are more likely to develop injection-related adverse events. However, patients with skin of color are more likely to develop post-inflammatory hyperpigmentation and/or hypertrophic scar/keloid formation following injection procedures. Patients with specific ethnic characteristics, e.g. Asian population, should be informed of a higher risk of tissue reactions, e.g. itching, swelling, erythema, inflammation.

10. Assembly of needle to syringe





. For optimal use of **INOVOSENSE SHINE Plus**, it is important that the needle is properly connected to the syringe. See diagrams 1, 2, 3 and 4.

1. Firmly hold the glass cylinder of the syringe and the Luer-lock adaptor between the thumb and the index fingers.
2. Grasp the protective cap with the other hand and break the seal.
3. Push & twist the needle on the syringe until a resistance is felt. Do not over-tight. Overtightening of the needle may lead to the Luer-lock moving and dislodging from the syringe.

4. Keep holding the Luer-lock and remove the sheath from the needle.

11. Storage

Store between 2°C and 25°C. Protect from light freezing. Avoid mechanical shocks.

Marketing Authorization Holder:

MEDISA ARA GOSTAR IRANIAN

Number 21, Jooybar Alley, Shad St,

Molla Sadra St, Tehran, Iran

Manufacturer:

ARYA TAJHIZ IBN SINA

Number 21, 20 St, Kaveh Industrial City, Saveh, Iran

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